

## Preparation of a copper plate sample for scanning electron microscopy

*Imaging results of a pure copper plate sample prepared by broad-beam ion milling for scanning electron microscopy (SEM).*

### Mechanical preparation

A sample of pure copper plate (10 x 10 x 1.5 mm) was cut and polished mechanically; diamond lapping films were used to produce a 1  $\mu\text{m}$  final polish. The sample was observed under SEM at 5 kV, SED (ET) at 1000x magnification (top image).

### Ion beam preparation

The sample was prepared using broad-beam ion milling with the Fischione Instruments Model 1060 SEM Mill.


The ion milling parameters used were:

Parameter	Setting
Voltage	4 kV
Beam angle	5°
Stage rotation	Rocking
Milling time	40 minutes

### Results and analysis

The sample was observed under SEM at 1 kV, SED (ET), 1000x magnification (middle image) and 5 kV, BSD, 1000x magnification (bottom image).

A significant amount of subsurface deformation may remain in the copper following mechanical polishing. Ion milling at a low incident angle with sample rotation removes this residual damage and produces a uniform



### COPPER SAMPLE OBSERVED UNDER SEM

*Top: 5 kV, SED (ET), 1000x. Middle: 1 kV, SED (ET), 1000x. Bottom: 5 kV, BSD, 1000x.*

planar surface ready for EBSD. The optimal choice of incident angle can produce electron channeling contrast at the same time, thus highlighting the grain and twin boundaries.