

Model	MSIL3	MSIL4	MSIL5	MSIL6	MSIL7
Noise reduction	12dB	15dB	15dB	12dB	12dB
Designed for	Small rotary pump	Single rotary pump	Large single rotary pump	Two rotary pumps	Water chiller
Compatible with	Pfeiffer UNO 5	Edwards RV8 / RV12 / RV18	Edwards E2M28	2 x Edwards RV12	V.d. Heyden 004
(similar to)	Pfeiffer DUO 5	Edwards XDS10 / nXDS10	Adixen 2033SD	2 x Adixen 2010L	
	Edwards RV5	Adixen/Varian/Agilent 2010L	Pfeiffer DUO 20M		
	Alcatel BB2002				
	Edwards E2M1.5				
Fans	2	2	2	4	4
Voltage	230/115V	230V / 115V	230V / 115V	230V / 115V	230V / 115V
Power rating	Fused 2300W	Fused 2300W	Fused 2300W	Fused 2300W	Fused 2300W
Power connections	IEC sockets	IEC sockets	IEC sockets	IEC sockets	IEC sockets
Outside dimensions	350 x 400 x 600mm	500 x 500 x 800mm	500 x 500 x 900mm	500 x 500 x 800mm	800 x 900 x 850mm
Inside dimensions	250 x 350 x 500mm	400 x 450 x 700mm	400 x 450 x 800mm	440 x 450 x 740mm	530 x 580 x 740
Weight	16kg	20kg	23kg	20kg	35kg