



Cleanliness



CS G4

User-friendly
cleanliness sink for
accurate particle extraction

Specialized - User friendly - Ergonomic




Benefits:

- **Timesaving**
 - By blank tests
 - By extraction
- **Ergonomic**
 - Height adjustable for operator
 - Control rinsing flow by foot pedal
- **Cascade filter mount**
 - Cartridge for both D47 and D90 filters

Features:

- Touchscreen
- Flow control in bar or lpm
- Preset amount of rinsing fluid with automatic stop
- Reminder for service intervals
- Settings for either alkaline chem or degreaser solution




Model 10 : 1000 x 600mm (internal dim.)

Model 14 : 1400 x 800mm (internal dim.)

Model 18 : 1800 x 1000mm (internal dim.)

Complete provider of laboratory consumables



References

VOLVO

Construction Equipment



Donaldson®
Filtration Solutions


Hydra
Specma



BOSCH

rexroth

A Bosch Company



Cleanliness Cabinet CC-G4

Specially designed cabinet for particle extraction of industrial components according to ISO16232, ISO4406/4407, VDA19.1 and custom standards.

Developed with 50 years of expertise in metrology and 20+ years of focus in technical Cleanliness.

Why choose NYLI

- Design, development, assembly and quality check all made in our facilities in Sweden
- Installation, service, training and support by specialized engineers trained by NYLI
- NYLI have an accredited laboratory for particle testing working with customers cleanliness specifications daily in our own systems.
- CE and RoHS compatible

